

Accent *Prime*

Love Your Shape

Accentuating Your Results with Two Effective Applicators

Accent Prime™ enhances the effectiveness with two unique applicators, enabling powerful treatments at unprecedented speed.

Colibri | *Non-surgical skin lift and wrinkle contraction*

The only market solution offering dual ablative and thermal modes for the treatment of sagging eyelids and facial wrinkles

The Colibri applicator embodies **the highest peak power** solution on the market for non-surgical eyelid lift, periorbital & perioral fine lines & wrinkles, and Xanthelasma.

Targeting hard-to-reach and hard-to-treat anatomical areas such as periorbital and perioral, this solution can focus on saggy and lax excess skin, which in other cases may need surgical intervention and prolonged downtime.

Colibri utilizes Unipolar RF-based technology to produce Microplasma sparks, delivering energy via a specially shaped tip that creates controlled micro-perforations. Once delivered, the micro ablation damage leads to skin contraction, resulted in tighter and smoother skin. Colibri offers a highly effective alternative to more invasive treatments used today.



Before



After 3 Treatments


Courtesy of Dr. Ines Verner MD, Dermatologist, Verner Clinic.

“The Colibri applicator is proving to be a highly valuable asset and addition to my Accent Prime treatments. My patients benefit from very effective treatments and much smoother, younger looking skin – with significantly shorter downtime than other alternate solutions.”


Dr. Ines Verner MD, Dermatologist, Verner Clinic

SpeedSculpt | Face & Body Contouring for Extremely Small Areas *Super Small. Super Sensitive Sculpting.*

The SpeedSculpt applicator is based on a revolutionary Ultrasound technology designed to treat small areas as part of any complete treatment. A down-sized, focused version of Alma’s UltraSpeed and MiniSpeed applicators, it emits guided ultrasonic waves with a small-sized 44 mm diameter sonotrode. Through concentric profiles, for high-speed fat volume reduction and contouring it is uniquely designed for effective and safe treatments on small areas. This applicator is perfect for use on exceptionally small and hard to reach areas, as well as for patients with a narrow physical structure.



Before



After 2 Treatments

Courtesy of Alma's Clinical Department

Alma Lasers GmbH

Nordostpark 100-102
90411 Nuremberg, Germany
Tel. + 49 911 / 89 11 29-0
Fax + 49 911 / 89 11 29-99

Email: info@almalasers.com

www.almalasers.com

© 2019 Alma Lasers. All rights reserved.
Alma, Accent Prime, its logos are trademarks or registered trademarks of Alma Lasers. Product specifications are subject to change without notice.

PBAC28021901

 CONNECT WITH
ALMA LASERS

